


FORCEPS

Desjardin gall stone forceps

screw joint

SC-IN/001a small, 220mm long


SC-IN/001b large, 238mm long



Babcock intestinal forceps

straight, extra light


SC-GI/278 203mm long



Maingot intestinal forceps

straight, serrated tip, fenestrated grip

SC-IN/006 203mm long



Denis Browne intestinal forceps

straight, serrated tips, fenestrated jaws


SC-IN/007 178mm long



Lane stomach holding forceps

straight, serrated jaws, screw joint

SC-IN/008 203mm long




Lahey cholecystectomy forceps

curved, longitudinal serrations, heavier pattern

SC-IN/011 76mm blades, 197mm long

SC-IN/011a 225mm long




Lahey Sweet forceps


SC-IN/011b 200mm long

SC-IN/011c 225mm long

SC-IN/011d 230mm long




FORCEPS



Moynihan's cholecystectomy forceps

curved on flat, screw joint

SC-IN/012 215mm long




Lane intestinal forceps

110mm blades, longitudinal serrations

240mm long

SC-IN/013 straight

SC-IN/013b curved




Seton Pringle intestinal forceps

straight, longitudinal serrations

178mm long

SC-IN/014 fine blades 2.5mm wide



Kocher intestinal forceps

with light spring blades and longitudinal serrations

SC-IN/015a straight, 90mm blades, 210mm long

SC-IN/015b curved, 85mm blades, 210mm long

SC-IN/015c straight, 120mm blades, 250mm long

SC-IN/015d curved, 1150mm blades, 250mm long


SC-IN/015e straight, 150mm blades, 250mm long

SC-IN/015f curved, 145mm blades, 250mm long

FORCEPS AND CLAMPS

Lang Stevenson intestinal forceps

SC-IN/018 82mm jaw, keyed blades, 203mm long



Parker Kerr intestinal forceps/clamp

blades 89mm, 254mm long


SC-IN/017a curved

SC-IN/017c straight

Guard/shoe

SC-IN/017b curved, fenestrated inner or outer curve, 89mm

SC-IN/017d straight, 89mm




Heartman intestinal clamp

with diagonal serrations, 200mm long

SC-IN/015g straight

SC-IN/015h curved




Schoemaker-Ogilvie colectomy forceps

straight, hollow blades

SC-IN/020 171mm long

SC-IN/020a 203mm long



Payr intestinal clamp


crushing, lever action, longitudinal serrations,

SC-IN/003 extra small, 50mm blades, 150mm long


SC-IN/002a small, 65mm blades, 210mm long

SC-IN/002b medium, 100mm blades, 240mm long

SC-IN/002c large, 150mm blades graduated in inches, 360mm long



CLAMPS / MISCELLANEOUS




Doyen intestinal clamp

240mm long

SC-IN/004a straight

SC-IN/004b curved




Lane twin stomach clamp

longitudinal serrations

large, 150mm blades, 305mm long

SC-IN/009b straight

SC-IN/009a curved



Zachary Cope anastomosis clamp

SC-IN/010a set: 3 clamps, lever action closing device and Tommy bar, 233mm long

SC-IN/010b set of 3 clamps

SC-IN/010c Tommy bar for set, 103mm long



Bake's bile duct dilators

malleable stems, 203mm long

SC-IN/016a set of nine

SC-IN/016c 3mmø head

SC-IN/016d 4mmø head

SC-IN/016e 5mmø head

SC-IN/016f 6mmø head


SC-IN/016g 7mmø head

SC-IN/016h 8mmø head

SC-IN/016i 9mmø head

SC-IN/016j 10mmø head

SC-IN/016k 11mmø head



Yankauer suction tube

SC-IN/019 305mm long