

Symmetry Surgical® Skin Markers



Symmetry Surgical's single-use, sterile **Skin Markers** are designed to provide a permanent skin marking solution for intact skin that won't wash off while preparing a patient for a procedure.

- Sterile Skin Markers require no preoperative preparation
- Packaged with a flexible, plastic ruler that follows body contours for easy marking and reference photos
- Available in standard tip for broad clear lines as well as in micro tips for extra precision
- Packaged 36 pre-sterile Skin Markers per box



Gentian Violet Dye, Sterile, Disposable, Packaging with 6 inch Ruler

| CODE | DESCRIPTION | QTY |
|---------|-------------|--------|
| 43-1030 | Regular Tip | 36/Box |
| 43-1032 | Micro Tip | 36/Box |

Distributed in the UK by



surgicalholdings

7-9 Parkside Centre, Potters Way, Southend-on-Sea
Essex SS2 5SJ T 01702 602050
E info@surgicalholdings.co.uk



symmetry surgical®

Instruments | Sterilization Solutions | Retractors | Energy

3034 Owen Drive | Antioch, TN 37013 | 800.251.3000 (US) | + 41 (0) 52 557 59 00 | customerservice@symmetrysurgical.com | symmetrysurgical.com

©2013 Symmetry Surgical, Inc. All rights reserved. Cat No. SM0042