



**surgical**holdings

UK manufacturer and repairer of Surgical Instruments

## SCISSORS

**UK manufactured  
quality scissors  
for general and  
plastic surgery**



Temple Farm Industrial Estate, Southend-on-Sea, Essex SS2 5SJ  
TELEPHONE: 01702 602050 FAX: 01702 460006  
info@surgicalholdings.co.uk www.surgicalholdings.co.uk



### **DUO SURFACE FINISH**

Surgical Holdings have utilised their combined expertise and knowledge to introduce a new and exclusive DUO surface finish.

By reviewing research and understanding the requirements of their customers they have combined two surface finishes on to one instrument. The first maintains anti-glare under theatre lights and the second ensures the best possible cleaning efficacy.

The company's highly skilled instrument finishers ensure that prior to any final stage finishing, all instruments receive a high level mopping to remove any imperfections.

From there they can then apply:

**The wet bead blast finish**

this technique uses ultra-fine glass beads sprayed at the surface of the instrument to create a matt looking finish.

**The mirror polish finish**

this technique uses a series of mops and compounds to polish the surface up to a mirror finish.

The result is a premium looking instrument which stands out, maintains functionally under theatre lights but above all has been proven to provide benefits against protein adhesion and corrosion.

For more detailed information ask for our [\*\*Premium Repair Service document\*\*](#)



Introducing **NERO**  
a new ceramic-coated  
range of scissors with  
unbeatable durability,  
life cycle cost,  
anti-glare,  
cleaning efficacy,  
and the smoothest  
ever performance.



Externally tested  
to 5000 usage cycles  
with no deterioration in  
cutting performance

## NERO SMOOTH CUT SCISSORS


---

### Metzenbaum scissors

curved on flat

**SC-GI/407bk** 185mm long

**SC-GI/408bk** 147mm long



---

### Mayo scissors

**SC-GI/410abk** curved on flat, 160mm long

**SC-GI/409bk** straight, 140mm long




---

### McIndoe scissors

178mm long

**SC-GI/150bk** curved on flat




## SUPERTWIST LIGHTWEIGHT SCISSORS

---




Introducing **SUPERTWIST**  
A lightweight scissor with a clever twist.  
The ergonomic design  
delivers improved dexterity, and  
is aesthetically innovative.

---



**Iris scissors**  
super cut  
114mm long


**SC-OP/476-st** straight  
**SC-OP/476a-st** curved



**Strabismus scissors**  
super cut  
114mm long

**SC-GI/475-st** straight  
**SC-GI/475a-st** curved


---



**Stevens tenotomy scissors**  
super cut, curved/blunt


**SC-OP/424-st** 114mm long  
**SC-OP/424a-st** 125mm long  
**SC-OP/424gb-st** 135mm long

---



**Jameson scissor**  
supercut, curved,  
micro-serrated cutting edge, 153mm long

**SC-SH/4179a-st**



**Kilner scissors**  
rounded blades, 115mm long

**SC-PL/001b-st** curved, blunt points  
**SC-PL/001c-st** straight, sharp fine points  
**SC-PL/001f-st** curved, sharp points

## PREMIUM QUALITY STAINLESS STEEL SCISSORS


---

### "Careless" scissors

165mm long

**SC-GI/474**

**SC-GI/474a** supercut




### Thomas Strabismus scissors

114mm long

**SC-OP/030** straight


**SC-OP/029** curved



### Cottle bulldog scissors

blunt/blunt, 115mm long

**SC-ENT/202** 115mm long




### Joseph scissors


sharp/sharp, 150mm long

**SC-ENT/203** straight

**SC-ENT/204** curved



## PREMIUM QUALITY STAINLESS STEEL SCISSORS




### Metzenbaum scissors

straight

- SC-GI/405** 180mm long
- SC-GI/406** 147mm long
- SC-GI/406a** 200mm long
- SC-GI/406b** 230mm long
- SC-GI/406c** 250mm long
- SC-GI/406d** 280mm long

curved on flat

- SC-GI/407** 185mm long
- SC-GI/408** 147mm long
- SC-GI/408a** 200mm long
- SC-GI/408b** 230mm long
- SC-GI/408c** 250mm long
- SC-GI/408d** 280mm long




### Mayo scissors

straight

- SC-GI/409** 140mm long
- SC-GI/409a** 160mm long
- SC-GI/411** 190mm long
- SC-GI/411a** 230mm long

curved on flat


- SC-GI/410** 140mm long
- SC-GI/410a** 160mm long
- SC-GI/412** 190mm long
- SC-GI/412a** 230mm long



### McIndoe scissors

178mm long


- SC-GI/149** straight
- SC-GI/150** curved on flat
- SC-GI/150sc** supercut, curved



### McIndoe scissors

curved on flat, round tapering ends

- SC-PL/007** 190mm long



### Potts scissors

satin finish, 190mm long

- SC-THO/019** angled on flat 45°
- SC-THO/020a** angled to side 45°
- SC-THO/020b** acute angle to side




**TUNGSTEN CARBIDE HARD INSERT**

**Mayo Harrington scissors**  
supercut, tungsten carbide

230mm long  
straight **SC-GI/162**  
curved **SC-GI/413**


200mm long  
straight **SC-GI/162b**  
curved **SC-GI/413b**



**Mayo scissors**  
supercut, tungsten carbide

straight **SC-GI/163** 140mm long  
**SC-GI/163a** 160mm long  
**SC-GI/165** 190mm long  
**SC-GI/165a** 230mm long


curved **SC-GI/164** 140mm long  
**SC-GI/164a** 160mm long  
**SC-GI/166** 190mm long  
**SC-GI/166a** 230mm long



**Metzenbaum dissecting scissors**  
supercut, tungsten carbide


straight **SC-GI/190tc** 147mm long  
**SC-GI/188tc** 185mm long  
**SC-GI/188tca** 200mm long  
**SC-GI/188tcb** 230mm long  
**SC-GI/188tcc** 250mm long  
**SC-GI/188tcd** 280mm long

curved **SC-GI/191tc** 147mm  
**SC-GI/189tc** 185mm  
**SC-GI/186tca** 200mm  
**SC-GI/189tcb** 230mm  
**SC-GI/189tcc** 250mm  
**SC-GI/189tcd** 280mm




**Gorney supercut scissors**  
(Rhytidectomy scissor) tungsten carbide,  
one micro serrated and one razor sharp edge

**SC-PL/085a** 125mm long  
**SC-PL/085** 200mm long, straight  
**SC-PL/085b** 200mm long, curved




**TUNGSTEN CARBIDE HARD INSERT** with two gold bows for identification




**Mayo dressing scissors**  
tungsten carbide hard insert

- |                   |            |                   |            |
|-------------------|------------|-------------------|------------|
| straight          |            | curved            |            |
| <b>SC-GI/419</b>  | 140mm long | <b>SC-GI/420</b>  | 140mm long |
| <b>SC-GI/419a</b> | 160mm long | <b>SC-GI/420a</b> | 160mm long |
| <b>SC-GI/421</b>  | 190mm long | <b>SC-GI/422</b>  | 190mm long |
| <b>SC-GI/421a</b> | 230mm long | <b>SC-GI/422a</b> | 230mm long |




**Mayo-Harrington dressing scissors**  
tungsten carbide hard insert

- |                  |          |                   |          |
|------------------|----------|-------------------|----------|
| 230mm long       |          | 200mm long        |          |
| <b>SC-GI/417</b> | straight | <b>SC-GI/162a</b> | straight |
| <b>SC-GI/418</b> | curved   | <b>SC-GI/413b</b> | curved   |




**Metzenbaum scissors**  
tungsten carbide hard insert

- |                   |            |                   |            |
|-------------------|------------|-------------------|------------|
| straight          |            | curved            |            |
| <b>SC-GI/445</b>  | 145mm long | <b>SC-GI/446</b>  | 145mm long |
| <b>SC-GI/448</b>  | 185mm long | <b>SC-GI/449</b>  | 185mm long |
| <b>SC-GI/471</b>  | 200mm long | <b>SC-GI/449a</b> | 200mm long |
| <b>SC-GI/471a</b> | 230mm long | <b>SC-GI/449b</b> | 230mm long |
| <b>SC-GI/471b</b> | 250mm long | <b>SC-GI/449c</b> | 250mm long |
| <b>SC-GI/471c</b> | 280mm long | <b>SC-GI/449d</b> | 280mm long |



**Baby Metzenbaum scissors**  
tungsten carbide,  
hard insert , 145mm long

- |                   |          |
|-------------------|----------|
| <b>SC-GI/445a</b> | straight |
| <b>SC-GI/446a</b> | curved   |




**Nelson's lobectomy scissor**  
tungsten carbide, curved on flat,  
230mm long

**SC-THO/057**

with two black bows for identification **SUPER CUT OPERATING SCISSORS**

**Metzenbaum dissecting scissors**

super cut



straight (length)

**SC-GI/190** 147mm  
**SC-GI/188** 185mm  
**SC-GI/188a** 200mm

curved (length)

**SC-GI/188b** 230mm  
**SC-GI/188c** 250mm  
**SC-GI/188d** 280mm

**Stevens tenotomy scissors**

super cut

curved/blunt

**SC-OP/424** 114mm long  
**SC-OP/424a** 125mm long  
**SC-OP/424gb** 135mm long



**SC-OP/426** straight/blunt, 114mm long  
**SC-OP/427** straight/sharp, 114mm long  
**SC-OP/427c** straight/sharp, 180mm long

curved/sharp

**SC-OP/425** 114mm long  
**SC-OP/425a** 140mm long  
**SC-OP/425b** 160mm long

with cutout, 110mm long

**SC-OP/427a** straight/blunt  
**SC-OP/427b** curved/blunt



**Stevens tenotomy 'Suregrip' scissors**


super cut, curved, bigger bow for extra support, control, and precision in use. Smoother action and improved dexterity

**SC-OP/424-dex** 125mm long

**Iris scissors**

super cut  
 114mm long

**SC-OP/476** straight  
**SC-OP/476a** curved



**Strabismus scissors**

super cut  
 114mm long


**SC-GI/475** straight  
**SC-GI/475a** curved



**Jameson scissor**


supercut micro-serrated cutting edge,  
 153mm long

**SC-SH/4179a**



**SUPER CUT OPERATING SCISSORS** with two black bows for identification

---



**Mayo scissors**

super cut

straight

**SC-GI/166a** 140mm long

**SC-GI/166e** 160mm long

**SC-GI/166b** 170mm long


curved

**SC-GI/166c** 140mm long

**SC-GI/166f** 160mm long

**SC-GI/166d** 170mm long

---




**Diethrich circumflex artery scissors**

super cut, 175mm long

**SC-THO/077** 125° angled on side

---





**Potts/Diethrich scissors**

super cut, 190mm long

**SC-THO/078** 25° angled on side

**SC-THO/079** 45° angled on side

**SC-THO/080** 90° angled on side



**Boyd scissor**

supercut, slight curve

**SC-THO/047** 180mm long

---



**Kilner scissors**

supercut, with black chrome bows

curved, blunt points,  
flattened at distal end.

**SC-PL/001e** 150mm long

**SC-PL/001d** 125mm long

straight, fine, sharp points  
with rounded blades,

**SC-PL/001c** 115mm long

---


## PLASTIC AND ORAL SCISSORS

---

### Cartilage cutting scissors

fine points, serrated edges


**SC-PL/009** 190mm long



### Fickling (Cathcart) scissors

angled

**SC-PL/049** 180mm long




### Kilner scissors

rounded blades, 115mm long

**SC-PL/001a** straight, sharp points

**SC-PL/001b** curved, blunt points




### Aufricht dissecting scissors

heavy pattern, triangular blades,  
140mm long

**SC-PL/002a** straight

**SC-PL/002d** curved




### Aufricht dissecting scissors

light pattern, 150mm long

**SC-PL/002b** straight

**SC-PL/002c** curved



## INDEX

---

<b>A</b>		<b>M</b>	
Aufricht dissecting scissors	9	Mayo dressing scissors	6
<b>B</b>		Mayo Harrington scissors	5
Baby Metzenbaum scissors	6	Mayo scissors	1,4,5,8
Boyd scissor	8	Mayo-Harrington dressing scissors	6
<b>C</b>		McIndoe scissors	2
"Careless" scissors	3	Metzenbaum dissecting scissors	5, 7
Cartilage cutting scissors	9	Metzenbaum scissors	1,4,6
Cottle bulldog scissors	3	<b>N</b>	
<b>D</b>		Nelson's lobectomy scissor	6
Diethrich circumflex artery scissors	8	NERO SMOOTH CUT SCISSORS	1
<b>F</b>		<b>P</b>	
Fickling (Cathcart) scissors	9	Potts scissors	4
<b>G</b>		Potts/Diethrich scissors	8
Gorney supercut scissors	5	<b>S</b>	
<b>I</b>		Stevens tenotomy 'Suregrip" scissors	7
Iris scissors	2,7	Stevens tenotomy scissors	2,7
<b>J</b>		Strabismus scissors	2,7
Jameson scissors	2,7	<b>T</b>	
Joseph scissors	3	Thomas Strabismus scissors	3
<b>K</b>			
Kilner scissors	2,8,9		

## QUALITY ASSURANCE

---

Surgical Holdings is an ISO 9001 and ISO 13485 accredited company. Its comprehensive quality management system is in accordance with the requirements of the Medical Device Directive 93/42/EEC - Annex I & II, for temporary orthopaedic surgical implants - skeletal pins, wires, staples, guide wires - and all associated instrumentation.

The company has compliance to manufacture and CE mark all Class I instrumentation, Class IIa and IIb implants, and can provide certificates of conformity for all materials used. All Surgical Holdings products undergo rigorous testing and inspection to applicable harmonised standards before leaving the factory. Complete traceability throughout all aspects of the production process provides the customer with access to a fully documented quality control system.

All instrumentation carries a lifetime guarantee from the date of purchase on any fault found as a result of poor materials. This warranty is void, however, if repairs are carried out by anyone other than Surgical Holdings repair personnel, or if it is deemed to have been misused.

Contact **01702 602050**  
or visit our website [www.surgicalholdings.co.uk](http://www.surgicalholdings.co.uk)  
for current reprocessing instructions



Temple Farm Industrial Estate, Southend-on-Sea, Essex SS2 5SJ

TELEPHONE: 01702 602050    FAX: 01702 460006

[info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)    [www.surgicalholdings.co.uk](http://www.surgicalholdings.co.uk)