




Single Patient Use

Ensuring a resilient supply chain of consumables with our trusted partners in Spain



Proud to be a family-run
British company

Our Spanish manufacturing partners 



[surgicalholdings.co.uk](https://www.surgicalholdings.co.uk)

Contents

Cover for Surgical Light Handle_4

Vein Strippers_5

Return Electrodes_6

Defibrillation Electrodes with Cable_10

Pre-Connected Defibrillation Electrodes_14

ECG Electrodes_17

Colourless Ultrasound Gel_19

Sterile Colourless Ultrasound Gel_19

Straight PVC Thoracic Catheters_20

Angled PVC Thoracic Catheters_21

PVC Trocar Thoracic Catheters_22

Vascular Silicone Ties_23

Yankauer Suction Set_24

Yankauer Cannulas_25

Symbol Key_26

If you have any questions about the products in this brochure,
please contact your Surgical Holdings representative for
more information on +44 (0)1702 602050



WE ALL HAVE AN
ETHICAL DUTY
TO REFURBISH
OUR MEDICAL INSTRUMENTS

Our European manufacturing partners offer reliability, speed, and sustainability

Product resilience and availability has been an issue since COVID in many areas. This has caused availability issues, cost escalation and transit challenges from the East. NHS Supply Chain established their resilience team in 2022 and we saw this as a key opportunity to see how we could help out and support this.

We decided to work with Telic and VEC Medical as key partners, due to their manufacturing being in Barcelona and their resilience focussed, single-use consumables. A reliable supply chain is ensured through availability and shorter delivery times due to the manufacturing location and the quick ease of transport from Spain, by road or sea.



We're proud to partner with Telic and VEC Medical to offer the high quality consumables shown in this brochure



Cover for Surgical Light Handle

Designed to prevent surgeon contact with the lamp handle during the surgical procedure.

The cover is sterile and completely covers the lamp handle.

Available in two sizes.



Reference	Code	Size	Compatible with handle	Units/Pouch	Units/Box
LHC-01	TCFOA	95x120mm	Ø 20-40mm	1	100
LHC-03	TCFOC	95x150mm	Ø 20-40mm	1	100




Vein Strippers

Dormo-Strip vein strippers, made entirely of polyamide and with an average tensile strength of 25kg for conventional or invagination stripping.

Each Dormo-Strip consists of 2 cables 100cm long, 3 different size olives and 1 handle.

Individually packed (1 bag, 1 kit), in boxes of 12 units.



Reference	Code	Procedure	Cable length	Cable diameter	Introducer diameter
VE-022	TCS01	Conventional	100cm	1.6mm	3.25mm
VE-025	TCS11	Invagination	100cm	1.6mm	3.25mm



RADIOPAQUE

STERILE

EO **Class IIa sterile**



Return Electrodes

Unipolar

Adult Unipolar

Total area: 217cm²

Contact surface area: 136cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2125	TPC01	-	-	×	1	100
2125-5	TPC05	-	-	×	5	100
2125	TPCT1	-	-	✓	1	100
2125-C/00	TPFOA	4200/3m	Erbe	×	1	50
2125-C/10	TPF0D	4210/3m	Valleylab	×	1	50
2125-C/00/5	TPFOG	4200/5m	Erbe	×	1	50



Paediatric Unipolar

Total area: 145cm²

Contact surface area: 83cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2225	TPD01	-	-	×	1	100
2225-5	TPD05	-	-	×	5	100
2225-C/00	TPGOA	4200/3m	Erbe	×	1	50
2225-C/10	TPG0D	4210/3m	Valleylab	×	1	50



Neonatal Unipolar

Total area: 83cm²

Contact surface area: 32cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2425	TPE01	-	-	×	1	50
2425-C/00	TPH0A	4200/3m	Erbe	×	1	50
2425-C/10	TPH0D	4210/3m	Valleylab	×	1	50

Dual

Adult Dual

Total area: 217cm²

Contact surface area: 128cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2500	TP101	-	-	×	1	100
2500-5	TP105	-	-	×	5	100
2500	TP1T1	-	-	✓	1	100
2500-C/00	TPN0A	4200/3m	Erbe	×	1	50
2500-C/12	TPN0E	4212/3m	Valleylab REM	×	1	50
2500-C/12/5	TPN0H	4212/5m	Valleylab REM	×	1	50



Return Electrodes

Dual

Adult Dual Oblong Shape

Total area: 221cm²

Contact surface area: 128cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2510	TPL01	-	-	×	1	100
2510-5	TPL05	-	-	×	5	100
2510	TPLT1	-	-	✓	1	100
2510-C/00	TPQ0D	4200/3m	Erbe	×	1	50
2510-C/12	TPQ0E	4212/3m	Valleylab REM	×	1	50

Paediatric Dual

Total area: 145cm²

Contact surface area: 73cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2600	TPJ01	-	-	×	1	100
2600	TPJT1	-	-	✓	1	100
2600-C/12	TPO0E	4212/3m	Valleylab REM	×	1	50
2600-C/00	TPO0F	4200/3m	Erbe	×	1	50





Neonatal Dual

Total area: 83cm²

Contact surface area: 31cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2700	TPK01	-	-	×	1	50
2700	TPKT1	-	-	✓	1	50
2700-C/00	TPP0A	4200/3m	Erbe	×	1	50
2700-C/12	TPP0E	4212/3m	Valleylab REM	×	1	50
2700-C/12	TPPTE	4212/3m	Valleylab REM	✓	1	50

Universal Dual

Total area: 168cm²

Contact surface area: 105cm²



Reference	Code	Connector/Cable	Compatibility	Traceability label	U/Pouch	U/Box
2900	TPT01	-	-	×	1	100
2900-5	TPT05	-	-	×	5	100
2900	TPTT1	-	-	✓	1	100
2900-C/00	TPV0A	4200/3m	Erbe	×	1	50
2900-C/12	TPV0E	4212/3m	Valleylab REM	×	1	50
2900-C/12/5	TPVOH	4212/5m	Valleylab REM	×	1	50


Defibrillation Electrodes with Cable




Desfi-Dormo defibrillation electrodes with cable are compatible with a wide variety of brands and models of defibrillation equipment.

The versions for adult patients are multifunction and can perform defibrillation, synchronised cardioversion and ECG monitoring procedures, although the compatibility with these functions will depend on the characteristics of the defibrillation equipment they work with.

Check the tables for the suitability of the connector with the equipment with which the electrodes must work and, if in doubt, consult your Surgical Holdings representative.

- Our electrodes are manufactured under the EN IEC 60601-2-4 standard.
- Supplied in packs of two defibrillator pads, in boxes of 10 packs.
- Available for adult and paediatric patients.



Connector	Reference	Code	Brand	Model	Patient	Multifunction
	EDC-1015	TO006	Physio-Control/Medtronic Mindray Lecor Bexen	Lifepack12, Lifepack15, Lifepack20e, Lifepack1000 Beneheart D3 Cor-Res A6S Reanibex 200	Adult	✓
	EDC-P115	TO051	Physio-Control/Medtronic Mindray	Lifepack12, Lifepack15, Lifepack20e Beneheart D3	Paediatric	✗
	EDC-1020	TO007	Agilent Philips	Heartstart XL, Heartstart FR2+	Adult	✓
	EDC-P120	TO052		Heartstart XL	Paediatric	✗
	EDC-1025	TO008	Agilent	Heartstream XLT (cable M3507A)	Adult	✓
	EDC-P125	TO053			Paediatric	✗
	EDC-1030	TO010	Progetti General Electric	Rescue Life, Rescue Sam Responder AED	Adult	✓
	EDC-P130	TO054	Progetti	Rescue Life vers. Display 7"	Paediatric	✗
	EDC-1035	TO012	Zoll M&B	Serie E, Serie R, Serie M AED7000, MB7000	Adult	✓
	EDC-P135	TO055	Zoll	Serie E, Serie R, Serie M	Paediatric	✗
	EDC-1040	TO011	Nihon Kohden	Cardiolife Tec-8300, Adapter cable JC-765V	Adult	✓
	EDC-P140	TO056		Cardiolife Tec-8300, Adapter cable JC-865V	Paediatric	✗

Connector	Reference	Code	Brand	Model	Patient	Multifunction
	EDC-1045	TO013	Schiller	Defigard 4000, Fred Easy	Adult	✓
	EDC-P145	TO057			Paediatric	✗
	EDC-1050	TO014	CU Medical Systems	Paramedic, CU-ER1	Adult	✓
	EDC-1055	TO015	Nihon Kohden	Cardiolife TEC-7731K, adapter cable JC-755V Cardiolife TEC-5521K, ActiBiphasic Cardiolife TEC-5531K, adapter cable JC-855V	Adult	✓
	EDC-P155	TO059			Paediatric	✗
	EDC-1060	TO016	HP / Philips Agilent Philips	Heartstart XL, Heartstart FR2+	Adult	✓
	EDC-P160	TO060	HP / Philips	Heartstart XL	Paediatric	✗
	EDC-1070	TO022	Corpuls	Corpuls 3 (except Kpome Slim)	Adult	✓
	EDC-P170	TO063			Paediatric	✗
	EDC-1090	TO025	Zoll	X Series	Adult	✓
	EDC-P190	TO066			Paediatric	✗


Pre-Connected Defibrillation Electrodes


The pre-connected Desfi-Dormo defibrillation electrodes allow faster and more effective action in case of emergency, as they can be directly connected to the defibrillator equipment.









The versions for adult patients are multifunction and can perform defibrillation, synchronised cardioversion, and ECG monitoring procedures, although the compatibility of these functions will depend on the characteristics of the defibrillation equipment they work with.

They are compatible with a wide variety of brands and models of defibrillation equipment. Check the tables for the suitability of the connector with the equipment with which the electrodes must work and, if in doubt, consult our commercial team.

- Our electrodes are manufactured under the EN IEC 60601-2-4 standard.
- Supplied in packs of two defibrillator pads, in boxes of 10 packs.
- Available for adult and paediatric patients.



Connector	Reference	Code	Brand	Model	Patient	Multifunction
	EDC-2015	TO00A	Physio-Control/Medtronic Mindray Lecor Bexen	Lifepack12, Lifepack15, Lifepack20e, Lifepack1000 Beneheart D3 Cor-Res A6S Reanibex 200	Adult	✓
	EDC-P215	TO05A	Physio-Control/Medtronic Mindray	Lifepack12, Lifepack15, Lifepack20e Beneheart D3	Paediatric	✗
	EDC-2020	TO00B	Agilent Philips	Heartstart XL, Heartstart FR2+	Adult	✓
	EDC-P220	TO05B		Heartstart XL	Paediatric	✗
	EDC-2035	TO00E	Zoll M&B	Serie E, Serie R, Serie M AED7000, MB7000	Adult	✓
	EDC-P235	TO05E	Zoll	Serie E, Serie R, Serie M	Paediatric	✗
	EDC-2040	TO00F	Nihon Kohden	Cardiolife Tec-8300, adapter cable JC-765V	Adult	✓
	EDC-P240	TO05F		Cardiolife Tec-8300, adapter cable JC-865V	Paediatric	✗
	EDC-2045	TO00G	Schiller	Defigard 4000, Fred Easy	Adult	✓
	EDC-P245	TO05G			Paediatric	✗
	EDC-2055	TO00I	Nihon Kohden	Cardiolife TEC-7731K, adapter cable JC-755V Cardiolife TEC-5521K, ActiBiphasic	Adult	✓
	EDC-P255	TO05I		Cardiolife TEC-5531K, adapter cable JC-855V	Paediatric	✗

Connector	Reference	Code	Brand	Model	Patient	Multifunction
	EDC-2060	T000J	HP / Philips Agilent Philips	Heartstart XL, Heartstart FR2+	Adult	✓
	EDC-P260	T005J	HP / Philips	Heartstart XL	Paediatric	✗
	EDC-P265	T005K	Mediana	D500	Adult	✗
	EDC-2070	T000L	Corpuls	Corpuls 3 (except Kpome Slim)	Adult	✓
	EDC-P270	T005L		Corpuls 3 (not Slim)	Paediatric	✗
	EDC-P275	T005M	Mediana	A10M-OK	Adult	✗
	EDC-2080	T000N	Schiller	Defigard touch 7 Fred PA-1	Adult	✓
	EDC-P280	T005N			Paediatric	✗
	EDC-2085	T000P	Telefunken Heartreset	HR1 FA1	Adult	✓
	EDC-2090	T000Q	Zoll	X Series	Adult	✓
	EDC-P290	T005Q			Paediatric	✗

ECG Electrodes


We have a wide range of Dermo electrodes for cardiac monitoring, which allows you to choose the most appropriate electrode for each circumstance, no matter how demanding it may be. Electrodes for short monitoring, stress tests and resting tests, and for adult, paediatric and neonatal patients.

For resting tests

Dermo-Tab resting electrodes help optimise testing costs.

With acrylic gel and a safe, soft and comfortable adhesion, which conforms to body shape, and its practical tab, which easily connects to the clamps, they are the best option for the diagnosis of patients at rest.

Radiotransparent.



Reference	Code	Patient	Measurements	Gel	Units/Pouch	Units/Box
T-2226	TEL01	Adult and Paediatric	22.5x35mm	Paper/PET with C-Ag/AgCl coating	100	5000



Class I



ECG Electrodes

For short monitoring with foam backing



EKF-22KT



SX-30



SF-40



LF-40



LF-50



SM-36

Dormo electrodes for short monitoring allow excellent results in short-term diagnostic tests. Their foam backing guarantees precise adhesion, and the stud connection and the solid or semi-liquid hydrogel texture adapt to any need in ECG tests.

Reference	Code	Patient	Measurements	Gel	Units/Pouch	Units/Box
EKF-22KT	TECA1	Neonatal	22x22mm	Solid	6	300
SX-30	TEM01	Adult, Paediatric, Neonatal	Ø 30mm	Solid	50	1000
SF-40	TEA12	Adult and Paediatric	32x40mm	Solid	100	2500
LF-40	TEAS2	Adult and Paediatric	32x40mm	Semi-liquid	100	2000
LF-50	TEAS1	Adult	Ø 50mm	Semi-liquid	50	1000
SM-36*	TEA05	Adult and Paediatric	36x42mm	Solid	50	1000

*Radiotransparent/magnetic resonance



Colourless Ultrasound Gel

Transonic Clear colourless gel for ultrasound diagnosis and therapy. Salt free. Does not irritate the skin. Non greasy. Hydrosoluble. Does not damage the transducer. It spreads out easily and uniformly.



Reference	Code	Content	Units/Box
G-15	TOG15	250ml	25
G-15/05	TOG24	500ml	20
G-15/1	TOG16	1000ml	20
G-15/5	TOG12	5000ml Flexible ¹	4
G-15/05	TOG24	500ml Rigid ²	4

1 Contains cannula and empty refillable 250ml bottle.
2 Contains dispensing pump and empty refillable 250ml bottle holding 4 units.

Class I

Sterile Colourless Ultrasound Gel



Reference	Code	Content	Units/Box
G-15/E	TOG26	20ml	48


STERILE R Class I sterile

Straight PVC Thoracic Catheters

Used in the thoracic cavity for air or fluid draining.

- Radiopaque line throughout the length of the tube
- Open and atraumatic distal
- Hole next to distal end breaks the radiopaque line as an XR reference
- Metric introduction scale with marks every 2cm

Code	Patient	Measurements
70010 PVC CH10	Paediatric	L 230mm Øin 2.2 * Øex 3.3mm
70012 PVC CH12	Adult	L 500mm Øin 2.5 * Øex 4.0mm
70016 PVC CH16	Adult	L 500mm Øin 3.7 * Øex 5.3mm
70018 PVC CH18	Adult	L 500mm Øin 4.2 * Øex 6.0mm
70020 PVC CH20	Adult	L 500mm Øin 5.2 * Øex 6.7mm
70024 PVC CH24	Adult	L 500mm Øin 5.8 * Øex 7.9mm
70028 PVC CH28	Adult	L 500mm Øin 6.8 * Øex 9.3mm
70032 PVC CH32	Adult	L 500mm Øin 8.0 * Øex 10.6mm
70036 PVC CH36	Adult	L 500mm Øin 9.0 * Øex 11.9mm
70040 PVC CH40	Adult	L 500mm Øin 10 * Øex 13.3mm




Angled PVC Thoracic Catheters

Designed for the post surgery (heart or thoracic) drain retrogressive insertion.

- Radiopaque line throughout the length of the tube
- Open and atraumatic distal
- Hole next to distal end breaks the radiopaque line as an XR reference
- Angled at 90°
- Metric introduction scale with marks every 2cm
- Sterilised with ethylene oxide

Code	Patient	Measurements
70120 ANGLED PVC CH20	Adult	L 550mm Øin 5.2 * Øex 6.7mm
70124 ANGLED PVC CH24	Adult	L 550mm Øin 5.8 * Øex 7.9mm
70128 ANGLED PVC CH28	Adult	L 550mm Øin 6.8 * Øex 9.3mm
70132 ANGLED PVC CH32	Adult	L 550mm Øin 8.0 * Øex 10.6mm
70136 ANGLED PVC CH36	Adult	L 550mm Øin 9.0 * Øex 11.9mm




PVC Trocar Thoracic Catheters

Designed for the post surgery (heart or thoracic) drain retrogressive insertion.

- Phtalate free Medical grade non-toxic PVC Catheters
- Aluminium and stainless steel
- Universal connection terminal bulb
- Metric introduction scale with marks every 2cm
- Smooth insertion distal
- Sterilised with ethylene oxide

Code	Patient	Measurements
70508 PVC CH08	Paediatric	L 230mm Øin 1.8 * Øex 2.7mm
70510 PVC CH10	Paediatric	L 230mm Øin 2.2 * Øex 3.3mm
70512 PVC CH12	Paediatric	L 245mm Øin 2.7 * Øex 4.0mm
70516 PVC CH16	Paediatric	L 300mm Øin 3.7 * Øex 5.3mm
70520 PVC CH20	Adult	L 400mm Øin 5.3 * Øex 6.7mm
70524 PVC CH24	Adult	L 400mm Øin 5.8 * Øex 7.9mm
70528 PVC CH28	Adult	L 400mm Øin 6.8 * Øex 9.3mm
70532 PVC CH32	Adult	L 400mm Øin 8.0 * Øex 10.6mm
70536 PVC CH36	Adult	L 400mm Øin 9.0 * Øex 11.9mm



Vascular Silicone Ties

Vascular Silicone Ties to facilitate a safe identification, retraction, occlusion, of veins, arteries, nerves, ureters and tendons during surgery. Delivering appropriate retraction and tension. Oval profile so as to avoid tissue damaging. No shape memory effects or creases.

Colour code ties: Blue = Veins, Red = Arteries, Yellow = Ureters, White = Nerves.

- Medical grade silicone
- Atraumatic oval shape
- Sterilised with ethylene oxide

Code	Measurements	Colour
95001	Mini: 1x2mm, 450mm length	Red
95002	Mini: 1x2mm, 450mm length	Blue
95003	Mini: 1x2mm, 450mm length	Yellow
95004	Mini: 1x2mm, 450mm length	White
95005	Maxi: 1.5x2.5mm, 450mm length	Red
95006	Maxi: 1.5x2.5mm, 450mm length	Blue
95007	Maxi: 1.5x2.5mm, 450mm length	Yellow
95008	Maxi: 1.5x2.5mm, 450mm length	White
95009	Super Maxi: 5.3x1.7mm, 450mm length	Blue
95010	Micro: 0.7x0.5mm, 450mm length	White



Yankauer Suction Set

Yankauer disposable Suction set reduces cross contamination and blocking risks. Intended to remove blood and similar liquids during surgery.

- Made of semi-rigid, malleable, non-toxic transparent PVC, which makes seeing the draining easier
- 2 to 4 open distal holes and atraumatic on the sides
- Anatomic handle, which makes it easier to hold it (comfortable grip)
- Great connection allows perfect adaptability to vacuum or suction tubes
- Two end options, with or without suction control
- Universal proximal end. 3 metre tube
- Sterilised with ethylene oxide

Code	Suction control
90412 Yankauer Cannula Set CH12	None
90414 Yankauer Cannula Set CH14	None
90416 Yankauer Cannula Set CH16	None
90418 Yankauer Cannula Set CH18	None
90420 Yankauer Cannula Set CH20	None
90422 Yankauer Cannula Set CH22	None
90428 Yankauer Cannula Set CH28	None
90512 Yankauer Cannula Set CH12	With control
90514 Yankauer Cannula Set CH14	With control
90516 Yankauer Cannula Set CH16	With control
90518 Yankauer Cannula Set CH18	With control
90520 Yankauer Cannula Set CH20	With control
90522 Yankauer Cannula Set CH22	With control
90528 Yankauer Cannula Set CH28	With control



Yankauer Cannulas

Yankauer is a suction handle with the purpose of removing the blood and similar liquids during surgery. It avoids major blocking risks.

- Made of semi-rigid, malleable, non-toxic transparent PVC, which makes seeing the draining easier
- 2 to 4 open distal holes and atraumatic on the sides
- Anatomic handle, which makes it easier to hold it (comfortable grip)
- Great connection allows perfect adaptability to vacuum or suction tubes
- Two end options, with or without suction control
- Universal proximal end
- Sterilised with ethylene oxide

Code	Suction control
90012 Yankauer Cannula CH12	None
90014 Yankauer Cannula CH14	None
90016 Yankauer Cannula CH16	None
90018 Yankauer Cannula CH18	None
90020 Yankauer Cannula CH20	None
90022 Yankauer Cannula CH22	None
90028 Yankauer Cannula CH28	None
90112 Yankauer Cannula CH12	With control
90114 Yankauer Cannula CH14	With control
90116 Yankauer Cannula CH16	With control
90118 Yankauer Cannula CH18	With control
90120 Yankauer Cannula CH20	With control
90122 Yankauer Cannula CH22	With control
90128 Yankauer Cannula CH28	With control



Symbol Key



Suitable only for high energies, up to 1100mA



Suitable for conventional procedures only



Compatible for connection with electro-surgical units with REM function



Does not contain natural rubber latex



Do not re-use



Sterilised using ethylene oxide



Sterilised by gamma radiation

RADIOPAQUE

Radiopaque



Time limit for use

Manufacturer and Repairer of Surgical Instruments

SMART REPAIR

SURGIFIX
ONLINE PORTAL

Smart Repair is our unique repair service for surgical instruments, offering the promise of long-term value for money in choosing the sustainable option, and a product replacement if we are unable to return the instrument to a 'as new' condition*.

In choosing our innovative Smart Repair service, customers benefit from dashboard access to our Surgifix Online Portal, enabling you to book, manage and monitor repairs online with your own assigned log in. Instantly communicate directly with the repair team to stay up-to-date with the live progress of repairs and the option to view historical scope repair reports.











Promoting education; sharing knowledge and offering training opportunities



*For the price of the repair. Applies to the 50 most used and repaired products, correct at the time of going to press.

Suppliers of

27

-  Surgical Holdings Instruments
-  Symmetry Surgical
-  Hermann Laparoscopic
-  Bissinger Bipolar
-  Light Leads
-  Din Baskets & Containers
-  Sterile Packed K-Wires
-  Custom Made Instruments

Refurbishment & Repairs of

-  Surgical Instruments
-  Laparoscopic Instrumentation
-  Orthopaedic Power Tools
-  Electrosurgery Devices
-  Rigid Endoscopes
-  Full Tray Refurbishment



Questions about any products in this brochure?
Contact your Surgical Holdings representative



surgicalholdings

Surgical Holdings Parkside Centre, Temple Farm Ind. Est., Southend-on-Sea, Essex, SS2 5SJ
E. info@surgicalholdings.co.uk T. +44 (0)1702 602050 F. +44 (0)1702 460006
surgicalholdings.co.uk



Version 1



Based in
the City of
Southend