


Single Use Vein Stripper

Built completely of polyamide with an extractor medium tensile strength 25kg. Design and finishing to facilitate introduction and sliding. The introducer diameter is 3.25mm and the cable diameter 1.5mm.

Each set consists of:

- 2 cables: 100cm length
- 3 olives: 9.5mm, 12.8mm and 15.4mm diameter
- 1 handle



STERILE EO

RADIOPAQUE



LATEX-FREE
KINDER TO STAFF
AND PATIENTS



PVC-FREE
KINDER TO THE
ENVIRONMENT

 Proud to be a family-run British company

E. info@surgicalholdings.co.uk

T. +44 (0)1702 602050

surgicalholdings.co.uk



Promoting education;
sharing knowledge and
offering training opportunities