

Cover for surgical light handle.

BLAYCO

> Product description

Cover for surgical light handle is intended to prevent the surgeon contacting intended or accidentally with the handle of the lamp during a surgical procedure.



> Product characteristics

- Prevents the surgeon from directly contact the surgical light handle during the surgical procedure.
- Suitable for all models of surgical lights featuring a handle with diameter within 20mm and 40mm.
- Latex Free.

> Technical specifications

> Materials

- Fixation discs
- Polyethylene homopolymer. (LDPE)
- Pouch:
Polyethylene.
- Packaging:scd
 - Inner microperforated PE pouch
 - Medical grade cellulose
 - Film Polyamide/Polypropylene

> Main dimensional properties

- External disc diameter: 11,8cm
- Pouch dimensions:
 - LHC-01: 9,5 x 12cm
 - LHC-03: 9,5 x 15cm
- Accepted diameters for light handle: 20 – 40mm

> Sterilization and expiry date

- Sterile product by ethylene oxide.
- Shelf life: 3 years upon the manufacturing date.

> Storage and maintenance

- Keep protected from direct sunlight / the outside.
- Keep in a dry place.
- Single use product. No maintenance required.

> Packaging

Ref.	Code	Presentation
LHC-01	TCF0A	Individual pouch. Service box of 100 units.
LHC-03	TCF0C	Individual pouch. Service box of 100 units.

Date: 2023 February

Specification N°: ETP/250/14

> Regulatory information

TELIC, S.A.U. guarantees that this product is in conformity with **Regulation 2017/745** and that it has been manufactured following the directives of the Quality Assurance System certified as **ISO 13485**.

This product is classified as:

- **Class I sterile** product according to Annex VIII of the Regulation (UE) 2017/745 rule 1.
- GMDN Code: 44977 – Operating/examination/treatment light handle cover, single-use
- EMDN code: V9099 (Various devices not included in other classes – Other).