

Medical devices and consumables for  
hospitals and healthcare centres



**blayco**<sup>®</sup> **DORMO**<sup>®</sup> TRANSONIC **+B.O**<sup>®</sup>



Distributed by: Surgical Holdings,  
Parkside Centre, Temple Farm Ind. Est.,  
Southend-on-Sea, Essex, SS2 5SJ, UK

Our experience of more than 40 years in the health sector, our commitment to our clients and our international presence in more than 100 countries, make us a trusted ally capable of offering the essential quality guarantee in the field of health.

We have a wide range of solutions adapted to the needs of medical professionals and we manufacture our products under strict quality controls.



---

Certified company



# INDEX

## **SURGERY** **4**

<b>Electrosurgery</b>	<b>5</b>
Electrosurgical plates	5
Electrosurgical pencils and accesories	11
<b>Cover for surgical light handle</b>	<b>15</b>
<b>Skin markers</b>	<b>16</b>
<b>Vein strippers</b>	<b>17</b>
<b>Protective pads</b>	<b>18</b>

## **CARDIOLOGY** **19**

Defibrillation electrodes	20
ECG electrodes	27

## **OTHER SPECIALITIES** **32**

Neurophysiology electrodes	33
Dermatological punches and currettes	38

## **NURSING** **39**

Ultrasound gels	40
Conductive gel	42
Lubricating gel	42
Nasal holder	43
Bite block	43
Otoscope speculums	44

## **PATIENT WELLNESS** **45**

Pre-gelled electrodes	46
Reusable pack for cold/hot	48
Paraffin	49
Protection and disinfection	49

# Solutions for Surgery

---



We offer a wide range of single-use solutions for the operating room environment, which guarantee maximum efficiency for the healthcare professional, as well as optimal and respectful handling of the patient.

---



# Electrosurgery

## Electrosurgical plates

### Blayco

We were the first manufacturers of electrosurgical plates in Spain, and to this day, we continue to be the only ones. Our long experience means that Blayco electrosurgical plates are designed and sized to offer all safety guarantees during use.

We have pregelled electrosurgical plates for adult, paediatric and neonatal patients, as well as models for universal use, unipolar or bipolar, with cable compatible with Erbe, Valleylab and REM Valleylab devices and without cable for being used with reusable cables.

There is a Blayco electrosurgical plate for every need that may arise in the surgical environment, no matter how demanding it may be.



**blayco**

24 Modelos de placas electrosurgidas:

1. Placa unipolar para adultos
2. Placa unipolar para pediatría
3. Placa unipolar para neonatos
4. Placa bipolar para adultos
5. Placa bipolar para pediatría
6. Placa bipolar para neonatos
7. Placa unipolar para adultos con cable compatible con Erbe
8. Placa unipolar para adultos con cable compatible con Valleylab
9. Placa unipolar para adultos con cable compatible con REM
10. Placa unipolar para adultos sin cable
11. Placa unipolar para pediatría con cable compatible con Erbe
12. Placa unipolar para pediatría con cable compatible con Valleylab
13. Placa unipolar para pediatría con cable compatible con REM
14. Placa unipolar para pediatría sin cable
15. Placa unipolar para neonatos con cable compatible con Erbe
16. Placa unipolar para neonatos con cable compatible con Valleylab
17. Placa unipolar para neonatos con cable compatible con REM
18. Placa unipolar para neonatos sin cable
19. Placa bipolar para adultos con cable compatible con Erbe
20. Placa bipolar para adultos con cable compatible con Valleylab
21. Placa bipolar para adultos con cable compatible con REM
22. Placa bipolar para adultos sin cable
23. Placa bipolar para pediatría con cable compatible con Erbe
24. Placa bipolar para pediatría con cable compatible con Valleylab

1. Placa unipolar para adultos  
2. Placa unipolar para pediatría  
3. Placa unipolar para neonatos  
4. Placa bipolar para adultos

**THICKNESS HYDROGEL 0,69mm**


**>15Kg**  
ADULT

**5-10Kg**  
PEDIATRIC

**50W MAX. <5Kg**  
NEONATAL


**>5Kg**  
UNIVERSAL

## Adult unipolar




Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2125	TPC01	-	-	✗	1	100
2125-5	TPC05	-	-	✗	5	100
2125	TPCT1	-	-	✓	1	100
2125-C/00	TPF0A	4200 / 3 m	Erbe	✗	1	50
2125-C/10	TPF0D	4210 / 3 m	Valleylab	✗	1	50
2125-C/00/5	TPF0G	4200 / 5 m	Erbe	✗	1	50

## Paediatric unipolar




Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2225	TPD01	-	-	✗	1	100
2225-5	TPD05	-	-	✗	5	100
2225-C/00	TPG0A	4200 / 3 m	Erbe	✗	1	50
2225-C/10	TPG0D	4210 / 3 m	Valleylab	✗	1	50

## Neonatal unipolar




83 cm<sup>2</sup> 32 cm<sup>2</sup>



Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2425	TPE01	-		✗	1	50
2425-C/00	TPH0A	4200 / 3 m	Erbe	✗	1	50
2425-C/10	TPH0D	4210 / 3 m	Valleylab	✗	1	50

## Adult dual




217 cm<sup>2</sup> 128 cm<sup>2</sup>



Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2500	TP101	-		✗	1	100
2500-5	TP105	-		✗	5	100
2500	TP1T1	-		✓	1	100
2500-C/00	TPN0A	4200 / 3 m	Erbe	✗	1	50
2500-C/12	TPN0E	4212 / 3 m	Valleylab REM	✗	1	50
2500-C/12/5	TPN0H	4212 / 5 m	Valleylab REM	✗	1	50



### Adult dual oblong shape



Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2510	TPL01	-		✗	1	100
2510-5	TPL05	-		✗	5	100
2510	TPLT1	-		✓	1	100
2510-C/00	TPQ0D	4200 / 3 m	Erbe	✗	1	50
2510-C/12	TPQ0E	4212 / 3 m	Valleylab REM	✗	1	50

### Paediatric dual




Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2600	TPJ01	-		✗	1	100
2600	TPJT1	-		✓	1	100
2600-C/12	TPO0E	4212 / 3 m	Valleylab REM	✗	1	50
2600-C/00	TPO0F	4200 / 3 m	Erbe	✗	1	50





## Neonatal dual



83 cm<sup>2</sup> 31 cm<sup>2</sup>



Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2700	TPK01	-		✗		
2700	TPKT1	-		✓		
2700-C/00	TPP0A	4200 / 3 m	Erbe	✗	1	50
2700-C/12	TPP0E	4212 / 3 m	Valleylab REM	✗		
2700-C/12	TPPTE	4212 / 3 m	Valleylab REM	✓		

## Universal dual



168 cm<sup>2</sup> 105 cm<sup>2</sup>






Reference	Code	Connector/Cable	Compatibility	Traceability label	Units/Pouch	Units/Box
2900	TPT01	-		✗	1	100
2900-5	TPT05	-		✗	5	100
2900	TPTT1	-		✓	1	100
2900-C/00	TPV0A	4200 / 3 m	Erbe	✗	1	50
2900-C/12	TPV0E	4212 / 3 m	Valleylab REM	✗	1	50
2900-C/12/5	TPV0H	4212 / 5 m	Valleylab REM	✗	1	50



## Reusable clamp-cable for electrosurgical plates

### Blayco

As a complement to the electrosurgical plates without cable, the three Blayco reusable clamp cable models available ensure the use of the plate with the majority of electrosurgical devices on the market.

	Reference	Code	Length	Connector	Plate type	Units/Pouch	Units/Box
	4200	TCC00	3 m	Erbe	Unipolar / Dual		
	4200	TCCA1	5 m				
	4210	TCC10	3 m	Valleylab	Unipolar	1	50
	4212	TCC12	3 m	Valleylab REM	Dual		
	4212	TCCA2	5 m				



## Electrosurgical pencils and accessories

### Single-use electrosurgical pencils

#### Blayco

Blayco single-use hand control electrosurgical pencils, easy to use, lightweight, with a comfortable anatomical profile, with a 3m or 5m cable.


Different options available, see the table below.

Individually packed in boxes of 50 units.



MB100

Reference	Code	Cable	Control	Presentation
MB100	TCM01	3 m	Fingertip	Blayco MB100 electrosurgical pencil + AB80 blade active electrode
MB100/5	TCM05	5 m	Fingertip	Blayco MB100 electrosurgical pencil + AB80 blade active electrode
MB200	TCM02	3 m	Fingertip	Blayco MB100 electrosurgical pencil + AB80 blade active electrode + AL-40 adhesive cleaning device for electrode



## Reusable electrosurgical pencil

### Blayco

Blayco MBR600 reusable electrosurgical pencil easy to use, lightweight, with comfortable anatomical profile, fingertip control, with 3 m cable and standard 3-pin connector.


Served with 70 mm blade active electrode, removable.



Reference	Code	Cable	Control	Presentation
MBR600F	TCR06F	3 m	Fingertip	Blayco MBR600 reusable electrosurgical pencil + 70mm blade active electrode



Clase IIb



## Cleaning device for electrode

### Blayco

Blayco AL-40 active electrode cleaning device, self-adhesive, sterile, radiopaque.

Reference	Code	Dimensions	Thickness	Units / Pouch	Units / Box
AL-40	TCG01	50 x 50 mm	6 mm	1	250



STERILE EO

RADIOPAQUE

Clase I estéril



## Active electrodes for electrosurgical pencil Blayco

Blayco active electrodes allow cutting and coagulation during electrosurgical procedures. They are compatible with all electrosurgical pencil models that have a standard 2.36mm connection.

Available in blade, needle and ball formats, tungsten or with non-stick coating to respond to any need in electrosurgery.



	Reference	Code	Length	Electrode	Units / Pouch	Units / Box	
	AB-50	TCE01	80 mm	Ball	1	50	
	AB-60	TCE02	70 mm	Needle			
	AB-70	TCE03	150 mm				
	AB-80	TCE04	70 mm	Blade			
	AB-90	TCE05	160 mm				
	ABC-45 HT	TCEB1HT	45 mm	Needle	5	500	TUNGSTEN
	ABC-55 HT	TCEB2HT	55 mm				
	ABC-65 HT	TCEB3HT	65 mm				
	ABT-50	TCET1	80 mm	Ball	1	50	NON-STICK COATED
	ABT-60	TCET2	70 mm	Needle			
	ABT-70	TCET3	150 mm				
	ABT-80	TCET4	70 mm	Blade			
	ABT-90	TCET5	160 mm				

## Cover for surgical light handle Blayco

Designed to prevent surgeon contact with the lamp handle during the surgical procedure.


The cover is sterile and completely covers the lamp handle.

Available in two sizes.

Reference	Code	Dimensions	Compatible with handle	Units / Pouch	Units / Box
LHC-01	TCF0A	95 x 120 mm	Ø 20 - 40 mm	1	100
LHC-03	TCF0C	95 x 150 mm			



## Skin marker Blayco



Blayco skin marker, available in two models: single tip 1 mm thick and double tip 0.5mm and 1mm thick.

Resistant and long-lasting violet ink. The ink basically contains gentian violet, a non-toxic, non-allergenic, and non-irritating component.

For use in general surgery, plastic surgery, orthopaedic surgery, neurosurgery, cardiovascular surgery, and radiation therapy.

Includes 15cm ruler.



Reference	Code	Dimensions	Ruler	Units / Pouch	Units / Box
RQ-01TO	TCH02TO	Regular tip 1 mm	✓	1	100
RQ-05TO	TCH05TO	Twin tip 0,5 mm + 1 mm	✓		




## Vein strippers

### Dormo-Strip

Dormo-Strip vein strippers, made entirely of polyamide and with an average tensile strength of 25Kg for conventional or invagination stripping.

Each Dormo-Strip consists of 2 cables 100cm long, 3 different size olives and 1 handle.

Individually packed (1 bag, 1 kit), in boxes of 12 units.



Reference	Code	Procedure	Cable length	Cable diameter	Introducer diameter
VE-022	TCS01	Conventional	100 cm	1,6 mm	3,25 mm
VE-025	TCS11	Invagination			

## Protective pad for surgical procedures

### Blayco-Pad

Blayco-Pad protective pad is a transpirable polyurethane foam pad coated with a medical grade acrylic adhesive mainly intended to prevent pressure sores.

It is adaptable to any contour and can be cut at will to cover any need.

Measurements: 300 x 14 x 200 mm.


Reference	Code	Units / Pouch	Units / Box
AC-3020	TCXA1	5	50





# Solutions for Cardiology

---



Our extensive range of ECG electrodes for monitoring, diagnosis and defibrillation for adult, paediatric and neonatal patients cover the full range of applications, and their different backing materials and hydrogel consistencies make them suitable for any patient and circumstance.

---

## Defibrillation electrodes with cable

### Desfi-Dormo

Desfi-Dormo defibrillation electrodes with cable are compatible with a wide variety of brands and models of defibrillation equipment.

The versions for adult patients are multifunction and can perform defibrillation, synchronized cardioversion and ECG monitoring procedures, although the compatibility with these functions will depend on the characteristics of the defibrillation equipment they work with.

Check the tables for the suitability of the connector with the equipment with which the electrodes must work and, if in doubt, consult our commercial team.

Our electrodes are manufactured under the EN IEC 60601-2-4 standard.

They are served in packs of two defibrillator pads, in boxes of 10 packs.

Available for adult and paediatric patients.



Connector	Brand	Model	Multifunction	Reference / Code Adult	Reference / Code Paediatric
	Physio-Control/ Medtronic	Lifepack12, Lifepack15, Lifepack20e Lifepack1000	✓	EDC-1015 TOO06	
	Mindray Lecor Bexen	Beneheart D3 Cor-Res A6S Reanibex 200			
	Physio-Control/ Medtronic	Lifepack12 Lifepack15 Lifepack20e	✗		EDC-P115 TOO51
	Agilent Philips	Heartstart XL Heartstart FR2+	✓	EDC-1020 TOO07	
	Agilent Philips	Heartstart XL	✗		EDC-P120 TOO52
	Agilent	Heartstream XLT (cable M3507A)	✓	EDC-1025 TOO08	
			✗		EDC-P125 TOO53
	Progetti	Rescue Life Rescue Sam	✓	EDC-1030 TOO10	
	General Electric	Responder AED			
	Progetti	Rescue Life vers. Display 7"	✗		EDC-P130 TOO54
	Zoll	Serie E, Serie R Serie M	✓	EDC-1035 TOO12	
	M&B	AED7000 MB7000			
	Zoll	Serie E Serie R Serie M	✗		EDC-P135 TOO55
	Nihon Kohden	Cardiolife Tec-8300 adapter cable JC-765V	✓	EDC-1040 TOO11	
		Cardiolife Tec-8300 adapter cable JC-865V	✗		EDC-P140 TOO56

Connector	Brand	Model	Multifunction	Reference / Code Adult	Reference / Code Paediatric
	Schiller	Defigard 4000 Fred Easy	✓	EDC-1045 TOO13	
			✗		EDC-P145 TOO57
	CU Medical Systems	Paramedic CU-ER1	✓	EDC-1050 TOO14	
	Nihon Kohden	Cardiolife TEC-7731K adapter cable JC-755V	✓	EDC-1055 TOO15	
		Cardiolife TEC-5521K ActiBiphasic	✗		EDC-P155 TOO59
	HP / Philips Agilent Philips	Heartstart XL Heartstart FR2+	✓	EDC-1060 TOO16	
	HP / Philips	Heartstart XL	✗		EDC-P160 TOO60
	Corpuls	Corpuls 3 (except Kpome Slim)	✓	EDC-1070 TOO22	
	Corpuls	Corpuls 3 (not Slim)	✗		EDC-P170 TOO63
	Zoll	X Series	✓	EDC-1090 TOO25	
			✗		EDC-P190 TOO66

# Pre-connected defibrillation electrodes

## Desfi-Dormo

The pre-connected Desfi-Dormo defibrillation electrodes allow faster and more effective action in case of emergency, as they can be directly connected to the defibrillator equipment.


The versions for adult patients are multifunction and can perform defibrillation, synchronized cardioversion, and ECG monitoring procedures, although the compatibility of these functions will depend on the characteristics of the defibrillation equipment they work with.

They are compatible with a wide variety of brands and models of defibrillation equipment. Check the tables for the suitability of the connector with the equipment with which the electrodes must work and, if in doubt, consult our commercial team.

Our electrodes are manufactured under the EN IEC 60601-2-4 standard.

They are served in packs of two defibrillator pads, in boxes of 10 packs.

Available for adult and paediatric patients.





Connector	Brand	Model	Multifunction	Reference / Code Adult	Reference / Code Paediatric
	Physio-Control/ Medtronic	Lifepack12, Lifepack15, Lifepack20e Lifepack1000	✓	EDC-2015 TOO0A	
	Mindray Lecor Bexen	Beneheart D3 Cor-Res A6S Reanibex 200			
	Physio-Control/ Medtronic	Lifepack12 Lifepack15 Lifepack20e	✗		EDC-P215 TOO5A
	Mindray	Beneheart D3			
	Agilent Philips	Heartstart XL Heartstart FR2+	✓	EDC-2020 TOO0B	
	Agilent Philips	Heartstart XL	✗		
	Zoll	Serie E, Serie R Serie M	✓	EDC-2035 TOO0E	
	M&B	AED7000 MB7000			
	Zoll	Serie E Serie R Serie M	✗		EDC-P235 TOO5E
	Nihon Kohden	Cardiolife Tec-8300 adapter cable JC-765V	✓	EDC-2040 TOO0F	
		Cardiolife Tec-8300 adapter cable JC-865V	✗		
	Schiller	Defigard 4000 Fred Easy	✓	EDC-2045 TOO0G	
	Nihon Kohden	Cardiolife TEC-7731K adapter cable JC-755V	✓	EDC-2055 TOO0I	
		Cardiolife TEC-5521K ActiBiphasic Cardiolife TEC-5531K adapter cable JC-855V	✗		

Connector	Brand	Model	Multifunction	Reference / Code Adult	Reference / Code Paediatric
	HP / Philips Agilent Philips	Heartstart XL Heartstart FR2+	✓	EDC-2060 TOO0J	
	HP / Philips	Heartstart XL	✗		EDC-P260 TOO5J
	Mediana	D500	✗		EDC-P265 TOO5K
	Corpuls	Corpuls 3 (except Kpome Slim)	✓	EDC-2070 TOO0L	
	Corpuls	Corpuls 3 (not Slim)	✗		EDC-P270 TOO5L
	Mediana	A10M-OK	✗		EDC-P275 TOO5M
 	Schiller	Defigard touch 7 Fred PA-1	✓	EDC-2080 TOO0N	
			✗		EDC-P280 TOO5N
	Telefunken/ Heartreset	HR1 FA1	✓	EDC-2085 TOO0P	
	Zoll	X Series	✓	EDC-2090 TOO0Q	
			✗		EDC-P290 TOO5Q

## Manual defibrillation electrode ED-1010 Desfi-Dormo

Pack of 2 Desfi-Dormo ED-1010 manual defibrillation and programmed cardioversion electrodes for adult patients.

Each package contains two white flexible foam electrodes with a 135 cm<sup>2</sup> acrylic-based hydrogel conductive area.

Our electrodes are manufactured under the EN IEC 60601-2-4 standard.

Reference	Code	Units / Pouch	Units / Box
ED-1010	TOO01	2-unit pack	50




## ECG electrodes

We have a wide range of Dormo electrodes for cardiac monitoring, which allows you to choose the most appropriate electrode for each circumstance, no matter how demanding it may be. Electrodes for short monitoring, stress tests and resting tests, and for adult, paediatric and neonatal patients.

### ECG electrodes for short monitoring - Foam backing Dormo

Dormo electrodes for short monitoring allow excellent results in short-term diagnostic tests. Their foam backing guarantees precise adhesion, and the stud connection and the solid or semi-liquid hydrogel texture adapt to any need in ECG tests.

	Reference	Code	Patient	Measurements	Gel	Units / Pouch	Units / Box
	EKF-22KT	TECA1	Neonatal	22 x 22 mm	Solid	6	300
	SX-30	TEM01	Adult Paediatric Neonatal	Ø 30 mm	Solid	50	1000
	SF-40	TEA12	Adult Paediatric	32 x 40 mm	Solid	100	2500
	LF-40	TEAS2	Adult Paediatric	32 x 40 mm	Semiliquid	100	2000
	LF-50	TEAS1	Adult	Ø 50 mm	Semiliquid	50	1000
	SM-36*	TEA05	Adult Paediatric	36 x 42 mm	Solid	50	1000
Radiotransparent / Magnetic resonance							



## ECG electrodes - Paper backing

### Dormo

The paper backing of the Dormo ECG electrodes for tests that, due to their duration, require offering the patient a greater degree of comfort and convenience, guarantees unique flexibility and adaptability, as well as great adhesion and breathability.

Their semi-liquid hydrogel texture guarantees excellent signal transmission in any condition.

Stud connection.



Reference	Code	Patient	Measurements (mm)	Gel	Units / Pouch	Units / Box
LP-50	TEALP	Adult	Ø 50 mm	Semiliquid	50	1000

## ECG electrodes for stress tests

### Dormo

Dormo electrodes with foam backing for stress tests are designed to guarantee their adherence to the patients' skin, ensuring signal transmission.

With semi-liquid gel and stud connection, they are the best option for short-term tests in the most demanding conditions.



Reference	Code	Patient	Measurements (mm)	Gel	Units / Pouch	Units / Box
LEH-36	TEDS1	Adult	36 x 50 mm	Semiliquid	50	1000



## ECG electrodes for resting tests

### Dormo-Tab

Dormo-Tab resting electrodes help optimize testing costs.

With acrylic gel and a safe, soft and comfortable adhesion, which conforms to body shape, and its practical tab, which easily connects to the clamps, they are the best option for the diagnosis of patients at rest.

Radiotransparent.



Reference	Code	Patient	Measurements (mm)	Gel	Units / Pouch	Units / Box
T-2226	TEL01	Adult Paediatric	22.5 x 35 mm	Paper/PET with C-Ag/AgCl coating	100	5000


## ECG electrodes for neonatal patients

### Dormo


Dormo ECG electrodes for neonatal and paediatric patients are manufactured under the strictest quality controls to guarantee optimal monitoring, thanks to the foam backing and the highest quality solid gel.

Reference	Code	Measurements (mm)	Cable	Connection	Units / Pouch	Units / Box
KFS-140	TEK26	Ø 22 mm	60 cm	Ø 1,5 mm Female safety plug	12	60
KFS-150	TEK27	Ø 22 mm	60 cm	Ø 4 mm Female connector	12	60

KFS-140



KFS-150



# Solutions for other specialties

---



Our range of Dormo neurophysiology electrodes and Blayco dermatological punches and curettes complete our offer of solutions for hospitals and healthcare centres.

---

## Neurophysiology electrodes

The line of Dormo neurophysiology products from Telic comprises high-quality solutions, which guarantee professionals comfortable handling and reliable results and, for patients, completely safe tests.

All electrodes are disposable, thus minimizing the risk of cross contamination and ensuring an optimal signal.

### Subdermal needles

#### Dormo

Dormo range of disposable subdermal needle electrodes offers single, twisted cable, and corkscrew needle electrodes.

Designed to minimize the risk of post-operative neurological effects in patients, the polishing of the needle tip is optimized to ensure easy and painless skin penetration.



### Single subdermal needles

Reference / Code	Measurements (mm)	Cable length	Presentation	Units / Case	Units / Box
TN140381M	14 x 0,38 mm (28G)	1 m			
TN1403815M	14 x 0,38 mm (28G)	1,5 m	4 units of each colour: red, green, black, white, yellow and blue.	24	1152
TN1904515M	19 x 0,45 mm (28G)	1,5 m			

Each case contains 24 individually pouched units.



## Twisted subdermal needles



Reference / Code	Measurements (mm)	Cable length	Presentation	Units / Case	Units / Box
TN130382M	13 x 0,38 mm (28G)	2 m	5 colours (green, blue, red, orange, yellow) combined with white or black	40 (20 pairs)	1200
TN1303825M	13 x 0,38 mm (28G)	2,5 m			
TN180382M	18 x 0,38 mm (28G)	2 m	6 colours (green, blue, red, orange, yellow and purple) combined with white or black	24 (12 pairs)	1152



## Corkscrew subdermal needles

Reference / Code	Cable length	Presentation	Units / Case	Units / Box
TNCS12M	1,2 m	4 units of each colour: red, green, black, white, yellow and blue.	24	1152
TNCS15M	1,5 m			

Each case contains 24 individually pouched units.








## Concentric needles

### Dormo

Dormo disposable concentric needles offer maximum safety and reliability during EMG tests, guaranteeing excellent signal quality, as well as patient comfort, reducing pain during use.



	Reference / Code	Measurements (mm)	Colour	Units / Case	Units / Box
	TNCC2503R	25 x 0,30 mm (30G)	Red		
	TNCC37045V	37 x 0,45 mm (26G)	Green	25	1200
	TNCC50045A	50 x 0,45 mm (26G)	Blue		



STERILE EO Clase IIa

## Cable for concentric needles

### Dormo

Reusable shielded cable for electromyography, for connecting the concentric needles to the equipment.

Cable length 1.25 m.

Reference / Code

TNCBLOC125M

Clase I



Compatible with most of the EMG equipment and disposables in the market



## Cup electrodes Dormo

Dormo cup electrodes are in the shape of small spoons that capture and emit the electrical signal from the patient's head.


Signal clarity, efficiency and ease of use are the main advantages of Dormo cup electrodes.

They are designed with flexible wires and optimized connector for easy headbox insertion and extraction.



Reference / Code	Measurements (mm)	Cable length	Presentation	Units / Case	Units / Box
TNCUP101M	Ø 10 mm	1 m	5 colours: blue, red, yellow, white and green.	25	5000
TNCUP1015M	Ø 10 mm	1,5 m			
TNCUP61M	Ø 6 mm	1 m			
TNCUP615M	Ø 6 mm	1,5 m			

Each case contains 25 individually pouched units.



Class I

## Surface electrodes

### Dormo

Dormo surface electrodes are designed to ensure patient comfort, are flexible and can be repositioned.

Its optimized connector makes it easy to insert and remove the head.

They are manufactured with an Ag/AgCl sensor and conductive hydrogel.



Reference / Code	Measurements (mm)	Cable length	Presentation	Units / Pouch	Units / Case	Units / Box
TN4D201M	Ø 20 mm	1 m				
TN4D202M	Ø 20 mm	2 m	Cables: black, green, blue and red	4	40	1200



## Neutral return electrode

### Dormo

The return electrode concentrates, collects and safely conducts out of the body, the current applied to the patient during neurophysiology procedures.



Reference / Code	Measurements (mm)	Cable length	Presentation	Units / Pouch	Units / Case	Units / Box
TN40501M	40 x 50 mm	1 m				
TN40502M	40 x 50 mm	2 m	Green cable	1	20	600



# Dermatological punches and curettes

## Sterile biopsy punch

Blayco

Blayco biopsy punches are sterile, with an ergonomically designed round plastic handle and a finely ground seamless, stainless-steel blade to ensure safe, clean and easy sample extraction.



Reference / Code	Diameter	Units / Case	Units / Box
BPS1.0	1 mm	10	500
BPS2.0	2 mm		
BPS3.0	3 mm		
BPS4.0	4 mm		
BPS5.0	5 mm		
BPS6.0	6 mm		
BPS7.0	7 mm		
BPS8.0	8 mm		



## Sterile dermal curette

Blayco

Blayco dermatological curettes have an excellent cutting capacity that makes working with them easier, compared to conventional products.

They have an easy grip handle and sharp side indicator.



Reference / Code	Diameter	Units / Case	Units / Box
DCS2.0	2 mm	10	500
DCS3.0	3 mm		
DCS4.0	4 mm		
DCS5.0	5 mm		
DCS7.0	7 mm		



STERILE R Clase IIa

# Solutions for Nursing

---



At Telic we have a wide variety of products to facilitate the work of the nursing professional and that are, at the same time, very convenient for the well-being of the hospitalized patient or in the process of diagnosis.

---

# Ultrasound gels

## Blue ultrasound gel

### Transonic

Transonic blue gel for ultrasound diagnosis and therapy.

Salt free. Does not irritate the skin. Non greasy. Hydrosoluble.

Does not damage the transducer.

It spreads out easily and uniformly.

Reference	Code	Content	Units / Box
G-15	TOG02	250 ml.	25
G-15/05	TOG25	500 ml.	20
G-15/1	TOG04	1000 ml.	20
G-15/5	TOG05	5000 ml Flexible* <sup>1</sup>	4
G-15/5RB	TOG07	5000 ml. Rigid* <sup>2</sup>	4

\*<sup>1</sup> Contains cannula and empty refillable 250 ml bottle.

\*<sup>2</sup> Contains dispensing pump and empty refillable 250 ml bottle purchasing 4 units.



Clase I

## Colourless ultrasound gel Transonic

Transonic Clear colourless gel for ultrasound diagnosis and therapy.

Salt free. Does not irritate the skin. Non greasy. Hydrosoluble.

Does not damage the transducer.

It spreads out easily and uniformly.



Reference	Code	Content	Units / Box
GC-15	TOG15	250 ml.	25
GC-15/05	TOG24	500 ml.	20
GC-15/1	TOG16	1000 ml.	20
GC-15/5	TOG12	5000 ml Flexible* <sup>1</sup>	4
GC-15/5RB	TOG08	5000 ml. Rigid* <sup>2</sup>	4

\*<sup>1</sup> Contains cannula and empty refillable 250 ml bottle.

\*<sup>2</sup> Contains dispensing pump and empty refillable 250 ml bottle purchasing 4 units.

Clase I

## Sterile colourless ultrasound gel Transonic

Reference	Code	Content	Units / Box
G-15/E	TOG26	20 ml	48

STERILE R Clase I estéril





## Conductive gel



Highly conductive gel.

For short and long ECG monitoring and electrostimulation.

Does not damage terminals or contacts.

Non skin-irritating. Water soluble.

Reference	Code	Content	Units / Box
G-10	TOG01	250 ml	25

Clase I

## Lubricating gel

Dormo

Dormo lubricating water-soluble gel, adequate for hospital and nursing processes that require non-specific lubrication.

Does not damage equipment or plastic, metal or textile instruments.

Non irritating.

Reference	Code	Content	Units / Box
G-20	TOGH1	100 ml	24

Clase I



## Nasal holder for gastric catheters

### Dormo-Nas

Nasal holder for gastric catheters for adult patients, easy to place and great fixation power. Prevents decubitus ulcers in the nasal tissues and eliminates the risk of kinking in the catheters.

The depth of catheter insertion can be adjusted as required.



Reference	Code	Compatible with catheters	Units / Pouch	Units / Box
7500	TADE1	14, 16, 18	1	100



## Bite block

### Mordedor-Mo Dormo

Rubber device that prevents possible damage to the endotracheal tube due to patient bites. Minimizes the risk of mouth ulcers. It should be washed, but it is not necessary to replace it during intubation.

Two sizes available: for adult and paediatric patients.

Length: 90 mm.



Reference	Code	Diameter	Compatible with catheters	Patient	Units / Pouch	Units / Box
7600	TOM01	Ø 19 mm	8, 9, 10	Adult	1	50
7650	TOM02	Ø 15 mm	6, 7, 8	Paediatric	1	50



## Disposable otoscope speculum

### Dormo-Spec

Disposable otoscope speculums made of polypropylene for adult and paediatric patients. Different models available depending on size, type of patient and compatible equipment.

Each case contains 250 units. Each box contains 1000 cases.



	Reference	Code	Diameter	Patient	Compatibility
	4010	TOK01	Ø 2.5 mm	Paediatric	<b>Heine</b> K-100 Beta 100
	4020	TOK02	Ø 4 mm	Adult	
	4030	TOK03	Ø 4 mm	Adult	<b>Riester</b> E-Scope RI-Mini Ri-Scope L1/L2 Pen-Scope
	4040	TOK04	Ø 2.5 mm	Paediatric	
	4050	TOK05	Ø 4 mm	Adult	
	4060	TOK06	Ø 2.5 mm	Paediatric	
	4070	TOK07	Ø 2.5 mm	Paediatric	<b>Heine</b> Beta 400 Beta 200 K 180 Mini 3000
	4080	TOK08	Ø 4 mm	Adult	
	4090	TOK09	Ø 2.5 mm	Paediatric	<b>Riester</b> Ri-Scope L3
	4095	TOK10	Ø 4 mm	Adult	
	4100	TOK11	Ø 4 mm	Adult	<b>Welch-Allyn</b>



CHECK COMPATIBILITY WITH UNIT



# Solutions for Patient Wellness

---



Our range of wellness solutions helps professionals make therapies more bearable for patients, minimizing the impact of the hospital environment and humanizing treatments.

---

## Pre-gelled electrodes for electrical stimulation Dormo-Tens

### DT Series

Dormo-TENS DT series reusable conductive silicone electrode for muscular stimulation, with reusable adhesive hydrogel, for single patient use. Long lasting.

Available RT-Series hydrogel spares, easy and practical replacement.

2 mm socket connection.



Reference	Code	Measurements (mm)	Units / Pouch	Units / Box
DT-50	TTEB1	50 x 50 mm	4	80
DT-100	TTED1	100 x 50 mm	2	40

### Hydrogel spares



Reference	Code	Measurements (mm)	Units / Pouch	Units / Box
RT-50	TTRB1	50 x 50 mm	20	80
RT-100	TTRD1	100 x 50 mm	10	40

## ST Series

Dormo-TENS ST series reusable soft and flexible non-woven electrode for muscular stimulation, with adhesive semiliquid hydrogel, single-patient use.

With preattached leadwire to fit a 2 mm female connection.



Reference	Code	Measurements (mm)	Units / Pouch	Units / Box
ST-30R	TTEC3	Ø 30 mm	4	120
ST-50	TTEC1	50 x 50 mm	4	120
ST-100	TTEC2	50 x 100 mm	4	60



## Reusable pack for cold/hot

### Dormo

Dormo reusable hot / cold bag, which can be used in both hot and cold therapies.

Cold is suitable for swellings, haematomas, contusions, sprains, muscular or joint injuries.

In hot mode it relieves muscle and joint pain and is suitable for low back pain, torticollis, rheumatism or hypothermia.

Available in two sizes.



Reference	Code	Measurements (mm)	Units / Box
FC-02	TTG00	130 x 180 mm	20
FC-01	TTG01	130 x 260 mm	20



## Cover for cold/hot pack

### Dormo

White plush polyester cover to protect the skin from direct contact with the cover and prevent burns.

With adjustable elastic band with velcro to facilitate the fastening in the different areas of application.



Code	Measurements (mm)	Units / Pouch	Units / Box
TTG10	130 x 260 mm	50	500

## Paraffin

Paraffin blocks of 2 kg for rehabilitation individually packaged.

Melting point of 46°C-50°C.

Code	Kg	Units / Box
PL1	1	10
PL2	2	5
PL5	5	5



 Clase I

## Protection and disinfection

### Dermoprotective liquid soap +BO

Very gentle liquid soap for frequent use. Contains moisturizing ingredients.  
pH 5.5 formula.

Reference / Code	Content	Units / Box
TLJ03	500 ml.	20
TLJ01	1000 ml.	20



## Surgical masks

### Blayco

Blayco 3 ply surgical mask, without glass fibres.


Hypoallergenic and with very low resistance to breathing.

Nose bar adaptable and high filtration capacity (>98%).













Perfect fitting.

Reference	Code	Units / Case	Units / Box
BM-01	TBB01	50	2000

Clase I



# Symbol key

	Plate total area (cm <sup>2</sup> )
	Plate contact area (cm <sup>2</sup> )
	Suitable only for high energies, up to 1100mA
	Suitable for conventional procedures only
	Compatible for connection with electrosurgical units with REM function
	Does not contain natural rubber latex
	Do not re-use
	Sterilized using ethylene oxide
	Sterilized by gamma radiation
	Not sterile
	Radiopaque
	Time limit for use





**Telic, SAU**  
**Headquarters and factory**  
Molí d'en Barri, 7  
08415 Bigues i Riells  
Barcelona – Spain  
T. +34 938 656 125  
telic@telic.es

**Telic, SAU**  
**Commercial branch and warehouse**  
Rieral de Can Barri, 9 bis  
08415 Bigues i Riells  
Barcelona – Spain  
T. +34 938 656 125  
telic@telic.es

**Telic France, SAS**  
**French branch**  
Rue Jean Baptiste Biot, 379  
66000 Perpignan  
France  
T. +33 468 927 885  
mail@telicfrance.es



Distributed by: Surgical Holdings,  
Parkside Centre, Temple Farm Ind. Est.,  
Southend-on-Sea, Essex, SS2 5SJ, UK

**[www.telic.group](http://www.telic.group)**

---